


## Application Note 2:

### Effect of Pre-treatment on Dynamic CO<sub>2</sub> Breakthrough Curves and CO<sub>2</sub> adsorption capacity



Dynamic CO<sub>2</sub> Breakthrough Curves of blank reactor (blue), 13X zeolite without (black) and with (red) heat-treatment. Gray dashed-line exhibits the gas switch from Ar to CO<sub>2</sub>. The red circles on each curve indicate the points where the first CO<sub>2</sub> signal was detected by the NDIR Sensor.





The time interval between the dashed line and the red circle, during which no CO<sub>2</sub> signal is detected, demonstrates the filling of the available pore volume of each porous sorbent material. This interval reflects the CO<sub>2</sub> adsorption capacity and efficiency of the different zeolite samples.

If no pre-treatment is applied, all the pores and internal cages remain occupied by ambient water and other gases. Therefore, no room for CO<sub>2</sub> capture.

The delay between the dashed line and the red circle can be attributed to a larger pore volume, which allows more CO<sub>2</sub> to be adsorbed before it becomes detectable by the NDIR Sensor, and a stronger adsorption affinity, whereby CO<sub>2</sub> molecules are more effectively captured by the zeolite active sites.

### Supporting Figure for water desorption during heat treatment.



(TG-MS results of 13X zeolite (a) % weight loss during heating under N<sub>2</sub> up to 400°C. (b) MS data recorded by Specson portable RGA 200 gas analyzer. Heating ramp is 10°C/min.)

One should always keep in mind the possible weight loss during heat treatment. In this specific example, 13X zeolite loses its 20% of initial weight. This is very important to calculate CO<sub>2</sub> capture capacity by unit mass.

Water desorption has started around 100°C, exhibited a peak maxima around 180°C, and continues around 20 min up to 370°C. Ease of outgassing depends on the sample nature.