



SPECSON  
INSTRUMENTS







# HIGH SHEAR MIXER



## Specson HSM-055

The Specson HSM-055 is designed to handle a wide range of mixing operations with a capacity range of up to 12 liters. Its 0.5 kW motor operates with a 220-volt single-phase (110 volts optional) and 50 Hz power supply, with a nominal maximum speed of 9000 rpm. The electric rise and fall stand allows for effortless raising and lowering of the mixing unit via the push-button control on the control unit.



# Specson HSM-055 Dimensions



# Specson HSM-055 Data Sheet

## Technical Specifications

Capacity	Up to 12L
Speed	Up to 9000 RPM
Motor Power	0.5kW (+optional 0.75kW)
Weight	11.2kg
Display	7 inch LCD Touch Screen
Body Parts	All wetted parts are 316 stainless steel



## Contact us for the customized solutions.

Specson Online Analiz ve Proje Deney Tasarımı Ltd Şti  
Mustafa Kemal Mahallesi, Cozone Yerleşkesi No:280, ODTÜ Teknokent

+90 507 922 11 02

infospecson@gmail.com