




SPECSON
INSTRUMENTS



PHOCAs



Specson PHOCAs

The Specson PHOCAs series consists of customized photocatalytic reaction units designed for advanced research applications. These integrated systems feature a magnetic stirrer, various types and numbers of LED lamps with different power ratings and emission ranges, a peristaltic pump unit for sampling, temperature monitoring and control, and a 7-inch capacitive touch screen for digital control of the entire process.

SPECSON PHOCAs



Contact us for the customized solutions.

Specson Online Analiz ve Proje Deney Tasarımı Ltd Şti
Mustafa Kemal Mahallesi, Cozone Yerleşkesi No:280, ODTÜ Teknokent

+90 507 922 11 02

infospecson@gmail.com