




SPECSON
INSTRUMENTS


Magnetic Stirrer



Specson IMS-8


Specson Instruments provides highly sensitive and accurate magnetic stirrer with 8-channel independent control capability. It comes with 7" capacitive touchscreen and aluminum case. User can easily set RPM and timer for each channel for their synthesis protocol. You can also contact our team for the customized solutions

Specson IMS-8



Technical Specifications

Weight	7.5 kg
Size	600x113 mm
Display Screen	7 inch LCD Capacitive Touch Screen
Timer	More than 100 hours for each channel
Stirring Rate	0-1000 Rpm
Frame	Aluminum 5754
Power	220 VAC



Contact us for the customized solutions.

Specson Online Analiz ve Proje Deney Tasarımı Ltd Şti
Mustafa Kemal Mahallesi, Cozone Yerleşkesi No:280, ODTÜ Teknokent
0 507 922 11 02